


Lab Skills

Practice: Pipetting Small Volumes




B3 Summer Science Camp
at Olympic High School
2016

Pipetter types

- Serological and micropipettes are used to accurately transfer small liquid volumes (micro-liter to milli-liter) accurately and precisely.
 - Continuously adjustable
 - Can be set to any transfer volume within its range, which is from 10% of the marked volume up to but not beyond the marked volume.

Micropipettors: 0.0001-1.0 ml



Serological: 1-25ml



Serological pipettes



Manual Pump



Pipetter



Electronic Pump



Meniscus Watch: eye level, lowest part – estimate volumes between lines

The MicroPipetter - 1

Make sure you know how to

- Select the proper micropipetter to transfer a specified volume of sample
 - volume needed must be inside the range of the device, the lowest volume is $1/10^{\text{th}}$ of the maximum volume, which is how the pipette is labeled.
- Set a specified volume on the pipette volume indicator using the volume adjustment knob
 - Volume is read on the side window
 - Adjustment knob is either the top plunger itself or just below the plunger
- Read the volume setting in correct units
 - For most of the micropipettors the unit is a microliter (ul), one millionth of a liter.
 - The 1ml micropipetter (this is 1000ul) shows tens of microliters (20 means 200ul or 0.200ml, with the final amount estimated with a division between 20 and 21 that is really between 200 and 210.


The MicroPipetter - 2

Make sure you know how to


- Select the correct tips and properly seat them (tap *lightly* holding the device vertically)
- To use: pipetter is held vertically, plunger is depressed, tip end is placed *a little* below the liquid surface, liquid is drawn up slowly, tip is withdrawn from the liquid, moved to target container, plunger is depressed slowly, tip is dragged up along side of container.
- Properly eject the tip into a waste container.

Parts of the Micropipettor

Parts of the Automatic Pipettor



Pipette Geography



Operating the Micropipette

Step 1: Select the correct pipetter and set the volume



Operating the Micropipette

Step 2: Read the volume

		
<p>(a): P-20 Model $6.86 \text{ mL} = 0.00686$ or $6.86 \times 10^{-3} \text{ mL}$</p>	<p>(b): P-200 Model $132.4 \text{ mL} = 0.1324$ or $1.324 \times 10^{-1} \text{ mL}$</p>	<p>(c): P-1000 Model $262 \text{ mL} = 0.262$ or $2.62 \times 10^{-1} \text{ mL}$</p>


Performance

Volume range	Inaccuracy (E%)		Imprecision (CV%)	
	Min. vol.	Max. vol.	Min. vol.	Max. vol.
Micropipettes				
0.1 - 2 µL	< +/- 6.0 %*	< +/- 2.0 %	< 5.0 %*	< 1.5 %
0.5 - 10 µL	< +/- 2.5 %**	< +/- 1.0 %	< 1.8 %**	< 0.5 %
1 - 10 µL	< +/- 2.5 %	< +/- 1.0 %	< 2.5 %	< 0.7 %
2 - 20 µL	< +/- 2.5 %	< +/- 1.0 %	< 1.7 %	< 0.5 %
5 - 50 µL	< +/- 1.5 %	< +/- 1.0 %	< 1.0 %	< 0.5 %
10 - 100 µL	< +/- 1.5 %	< +/- 0.8 %	< 1.0 %	< 0.2 %
20 - 200 µL	< +/- 1.5 %	< +/- 0.8 %	< 0.6 %	< 0.2 %
100 - 1000 µL	< +/- 1.5 %	< +/- 0.5 %	< 0.5 %	< 0.2 %

* At 0.5 µL. Indicative data at 0.2 µL: E < +/- 12 %, CV < 8 %

Operating the Micropipette

Step 3: Attach the disposable tip




Attaching the disposable tip

Tip boxes are labeled for the range of barrel they accommodate (there are some manufacturer-specific types so check if you are mixing and matching).



Center the end of the barrel in the tip and tap straight down, gently, twice.

Some tips include a filter barrier – mostly used for PCR reactions.



Operating the Micropipette

Step 4: Depress the Plunger to the First Stop



Depress the plunger first: don't put the tip in the liquid and then depress the plunger (the pressure is different).

Step 5: Immerse Tip in Sample



The tip must be below the surface of the liquid throughout the process of pulling up the liquid.

Ideally the pipette is held vertically.

To aspirate the sample into the tip, allow the pushbutton to return *slowly and smoothly* to the fully extended starting position.

NEVER LET THE PLUNGER SNAP UP!

Leave the tip in the solution for 1-2 sec (longer for a viscous solution), then slowly withdraw.

Operating the Micropipette

Step 6: Withdraw the tip from the solution

Withdraw the pipet from the receiving vessel carefully, touching or sliding the tip along the wall of the vessel.





Note – you should be wearing gloves – no bare hands!

Operating the Micropipette

Remove the tip from the sample liquid. There should not be liquid on the outside of the tip if you have used it correctly – generally we do NOT touch the tip with anything. Some reagent containers may not allow ideal practice – in that case wipe away any droplets on the outside of the tip with a lint-free tissue, such as KIMWIPES. Don't touch the tip opening or you will wick away some of the sample. Note on the left that the volume is not correct (air at bottom).

Step 7: Withdraw the Tip



Note – you should be wearing gloves – no bare hands!

Operating the Micropipette

Step 8: Dispensing the solution to a tube



(a) Start
Dispensing

(b) 1st Stop =
Dispense

(c) 2nd Stop =
Expel



- Touch the tip end to the bottom or side wall of the receiving vessel
- Depress the plunger **slowly** to the FIRST STOP.
- Pause for at least one second, longer for larger volumes or viscous liquids.
- Press the plunger to the SECOND STOP (the second point, of greater resistance, at the bottom of the stroke) to expel any residual liquid in the tip – an air bubble should force out the last drop of liquid.


Operating the Micropipette

Step 9: Use the ejector button to eject the tip into a waste container (use a Biohazard can if you are directed, otherwise an empty tip box works pretty well).




First Lab: Pipetting Skills

- A set of pipette pumps. Serological pipettes, micropipettors and tips
- Several capped solution containers (colored water – why?, then glycerol – why?)
- Parafilm squares and small weigh boats
- Mass balances (2 levels – why?)



Checking the volume

- If you know the density (mass per volume) of the liquid you are pipetting, **and** you have a calibrated mass balance,
- You can pipette a specified volume into a pre-weighed (tared) container and determine the mass of the volume added
 - Why use water to check your pipetting?



At 20°C, H₂O has density as shown:

$$1\text{ml} = 0.9985\text{ gm}$$

$$1\text{ul} = 0.00099\text{ gm} = 0.99\text{ug}$$

A very close approximation of 1ug to 1ul

Using the Balance



Pick the right balance: does it go to a small enough amount (you need micrograms for most of the lab)?


Close the glass doors: Air currents can change the reading on a very sensitive balance.

Make sure you know where the Tare (Zero) control is.

Using Parafilm or weigh boats

Parafilm for volumes < 1ml

Use a weigh boat for volumes > 1ml



Aqueous solutions bead up on Parafilm – you can leave the backing paper on, but pipette solution onto the waxy side **not** the paper side.

After weighing the volume, blot it up from the Parafilm with a Kimwipe, re-tare and repeat the pipetting.



Accuracy and Precision

- Accuracy means the closeness with which the dispensed volume approximates the volume set on the pipette
- The level of accuracy is specified as mean error, the average deviation of *replicate measurements* from what is expected from the volume you set.
- Precision is the "scatter" or variance of individual measurements obtained from the same volume setting.
- Precision can also be expressed as standard deviation (variance divided by the mean).

Accuracy and Precision (Continued)

- Device capabilities: relative *accuracies* are generally about 1% or less for micropipettors within range
 - These micropipettors have recently been calibrated.
- Precision error is less than 0.5 % except when transferring the smallest *recommended* volume for a given pipette model
 - Using the pipettes to transfer volumes which are below the recommended range will introduce larger errors

Lab Practice with Pipettes



- Practice setting a few volumes
- Practice reading the digits of set volumes
- Practice seating the tip, drawing up and dispensing samples of water and of a glycerol solution (glycerol is denser *and* more viscous).
- Get the "feel" of the 1st (set volume to pull up) and 2nd (blow-out for delivery) stops
- Practice using the pipette and record how well measurements match settings (indirectly, your skill).